

## IVUS-guided Aortic Branch Cannulation – in vitro protocol

- Each phantom should take 5 mins per user - Order – Jonah, Boat, Brett
- for 3 users - that's ~10-15 mins per phantom
- Each phantom should take 15 mins to map
- **Things to point out when participants arrive** - best to go above and let **fall in**, reference point on the wire, fast as possible (speed), verifying cannulation with fluoro, don't overdrive, warning to not move bed, show steerable catheter, 90 degree view rotation, authorship

**K6 (pva mapped before, give demo), K1, K12, K5, K3, K4 (pva), K8 (pva)**

### Room Setup:

- Ensure **Nikon camera is set to 25 fps**, long lens, place very up close on high stand [ ]
- ensure correct ML model loaded [ ]
- presave data directories on computer **for all 7 phantoms and load preop data for each!** [ ]
- Have deformation registration directory also open for troubleshooting [ ]
- Ensure yaml parameters look correct [ ]
- Turn pump frequency down to ~70bpm [ ]
- Use syringe to remove air from pulsatile pump [ ]
- ensure tubings submerged where possible [ ]
- Lock the table in place [ ]
- Setup camera looking at fluoro screen and surgeon + fluoro + IVUS [ ]
- Reference point on wire [ ]

### Mapping of phantom (for each phantom):

- Place phantom in tank, remove bubbles
- If not enough motion, change construction
- Position fluoro
- (Inkbit) 16 Fr sheath
- (PVA) Feed guidewire up, Clip IVUS catheter into pullback device
- Start recording with (1) **cameras** in room, (2) **screen** recording, (3) click '**start recording**' in aortascope, and (4) click **start pullback device** in that order

### Registration contingency plans (20 mins)

1. **Run registration again** – or try viterbi
2. **Replay data** – conf class 2 threshold, tsdf voxel size, ML model
3. **Registration tuning** – dbscan eps, downsampling CT mesh, allowable stretch factor, angle filter ransac, renal order swap injective mappings, bin / full selected for 4D reg
4. **IVUS mapping only (say nthng)** and record transform data (superimpose CT scan later)

### Side branch navigation:

0. Cover phantom with black cloth, plug in steerable catheter
0. Check if wire is still decent quality / exchange at end of each run
0. Be sure to save previous surgeon's **transform data** (if using IVUS TSDF or deforming phantom)
0. Turn on **pulsatile pump** – for video even if pullback device not used
2. **record with** (1) **cameras** in room, (2) **screen** recording (above all else)
3. **Start the clock - manipulated to be at 90 degrees** to the side branch vessel for clear side view
  - perform a **shot of fluoroscopy to verify** successful cannulation, press **Lap** function
4. Every time surgeon **fires fluoro = 1 cannulation attempt**

**Materials checklist:**

Remember to book relevant materials on lab google calendars

**Computational cart:**

- tape
- IVUS machine
- EM station
- EM emitter
- EM tracker stand (Anthony lab)
- PIM module
- USB splitter
- marker for drawing on fluoro screen
- basin to prevent water leakage
- pump-aorta large inlet connector
- black cloth
- second cart for monitor (for now)
- pullback device
- pullback device power chord
- pullback device USB cable
- pulsatile pump
- computer
- computer monitor facing out
- computer charger
- frame grabber
- optional: frame grabber for fluoro
- optional: additional stand for computer monitor
- mouse
- keyboard
- supply boxes
- dash monitor
- camera
- camera stand
- iPhone stand
- data recording sheets
- pens
- screwdriver
- scissors
- lighter
- beaker
- extra basin for dumping water
- wooden/plastic board as metal table alternative
- barbell weights
- guidewires consolidated in box sizes
- loop closure connectors
- PVA phantom boxes
- inkbit phantoms
- M6 wrench
- towel
- dryseal

- velcro (male and female)

**supplies:**

- printed protocol
- barbell weights
- IVUS + EM catheters
- steerable EM catheter
- steerable catheter
- computer monitor
- guiding catheters
- guidewires
- contrast agent
- puffing catheters
- surgical drapes
- syringes for contrast agent injection

**FEVAR specific:**

- ensure surgeons wearing scrubs and lead
- add nester coils if needed, apply whiteout, add nester coils (record video), resheath graft
- map and register the aneurysm
- (record video) deploy graft while recording
- map graft with guidewire – change vpC, sim device, refine, deeplumen params as needed
- call centreline\_caller, rename refine
- (record video) register graft with beads
- (record video) fenestrate graft, clicking highest to lowest orifices
- cannulation
- **FEVAR branch deployment AFTER ALL studies**